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INTENDED USE

The dental burs are attach to the dental High/Low Speed Dental handpiece for use in dental clinics in removing carious materials, reducing hard tooth structure, cavity preparation, finishing tooth preparations and restoration and polishing teeth.

APPLICATION

Cavity preparation	Root canal preparation	Root planing	Crown and bridge technique
Crown preparation	Oral surgery	Removal of old fillings	Crown removal
Working on fillings	Model casting technique	Acrylic technique	Implantology
			Spray-cooling

SAFETY RECOMMENDATIONS

1. Be sure that the bur is inserted fully and tightly into the handpiece.
2. Select suitable speed according to recommended speed table. (Table 3)
3. Never exceed the maximum permissible speed. Not observing the maximum speed will result in the increased risk of injury.
4. Excessive contact pressure has to be avoided because excessive contact pressure may lead to the increased risk and the damage.
5. Using Z-Diamond burs in different materials: 1) material of Zirconia, including all diamond grits of all types; 2) other materials, especially coarse grit and super coarse grit, are recommended that the users observe the following information before use:
 - 5.1. Care should be taken to provide adequate cooling and minimal pressure. (Please always use spray cooling (minimum 50 ml/min)).

- 5.2. In order to reach an optimal roughness depth after the usage of a coarse or super coarse diamond instrument the subsequent use of a finishing bur might be essential.
6. Improper use leads to increased risk and inferior results. Please make sure the application and all instructions for use.
7. Discard any damaged bur or worn working parts immediately because the damaged or worn bur will lead to the user applying higher contact pressure which results in the increased operating temperature. This may lead to the increased risk and injury.

Table 3 Recommended Speed Ranges

Shank (FG)		Shank (RA)	
Head diameter in 1/10 mm	Maximum (RPM)	Head diameter in 1/10 mm	Maximum (RPM)
008 – 010	160,000	008 – 023	100,000
012 – 014	↑		
016 – 018	300,000	Shank (HP)	
021 – 023	↓	Head diameter in 1/10 mm	Maximum (RPM)
024 – 027	↓	012 – 033	100,000
031 – 032	160,000	040 – 053	↓ 80,000
033 – 040	↓	060 – 070	↓ 35,000
042 – 050	80,000		

STERILIZATION REQUIREMENTS/PROCEDURES RECOMMENDATIONS

- The dental bur is intended to be sterilized with a certified type B steam sterilizer
- The sterilization parameters are

Cycle type	Exposure Temperature	At least exposure Time
dynamic air removal	134°C	4 minutes


- Washing and rinsing is the first step. Wash in warm soapy water or a compatible detergent with a ultrasonic cleaner in 5 minutes. Rinse the dental burs thoroughly in clean water.
- Put the dental burs into the sterilization pouch together with a certified chemical index (CI), and then seal the pouch.
- Following the sterilizer user instruction to sterilize the sterilization load (sealed pouches).
- When completing the sterilization, unpack the pouch prior to use and verify if the CI is properly indicated.
- If re-use without sterilization, the risk in cross-infection and/or contamination could harm the patient.
- It is proved by experiments that diamond burs can still retain a certain effective cutting performance after 3 consecutive cycles of use. One cycle is cleaning and sterilizing diamond bur after 5 minutes of continuous use.

Warning


- The device is to be sterilized prior to be used; without the sterilization could cause bio-contamination and/or cross-infection that could harm the patient seriously.
- The device is intended to be re-use up to times; over-time use could degrade the intended performance.

Ordering options

- Dx Order No.:** prefer to use Dentex number system, in order to prepare orders efficiently by computer.
 Ex. 1.1. C 206 (*C : Diamond grit & designation、*206 : Dx No.) (Fig.1)
 1.2. C Z206 (*C : Diamond grit & designation、*Z206 : Dx No.- Zirconia) (Fig.2)
- Reference No.:** refer to the common use codes of European brands.
 Ex. 2.1. C 836 012 (*C : Diamond grit & designation、*836 : Reference No、*012 : size \emptyset 1/10 mm) (Fig.1)
 2.2. C Z836 012 (*C : Diamond grit & designation、*Z836 : Reference No、*012 : size \emptyset 1/10 mm) (Fig.2)
- ISO No. :** accord to ISO 6360 Number Coding System and apply to the abbreviated form

(Fig.1)

(Fig.2)

806 314 110 534 012
 ① ②③ ④ ⑤

- ① : Material of the working part**
 - Diamond, galvanic metal bond
- ②、③ : Shank and overall length**
 - FG、19 mm
 - Coupling dimensions according to ISO 1797 (Table 1)
- ④ : Shape and design (Cylindrical, side-and end-cutting / Coarse grit, hard bond)**
 - Diamond grit design according to ISO 7711-3 (Table 2)
- ⑤ : Nominal size ISO 2157**
 - Largest diameter of the working part (\emptyset 1/10mm)
 - Ex. : 012=1.2 mm

DENTEX DENTAL INDUSTRIAL CORP.
 8F., No. 8, Lane 348, Sec. 2, Chung Shang Rd.,
 Chung-Ho Dist., New Taipei City 235, Taiwan.

2460

DENTEX DENTAL INDUSTRIAL CORP.

E-mail: dentex@ms77.hinet.net

Tel: +886-2-22424333

Fax: +886-2-22424345

Website: www.dentex.com.tw

SMALTO

E-mail: smalto@otenet.gr

Tel: (+30)2310229686

Website: www.pastelli.gr

SMALTO
 Vas. Irakliou 19, 54625, Thessaloniki, Greece.

Version:2020/05